



**SCOTTISH
COCHLEAR IMPLANT
PROGRAMME**

Connecting to a Games Console

How to use an assistive listening device to connect to a Games Console



SCOTTISH
COCHLEAR IMPLANT
PROGRAMME

Before we begin, make sure you have a microphone paired to your processors.



Roger Device



Mini Mic

There are three different ways to connect to your games console. You can:

- Stream game audio to your processors **ONLY**.
- Use a headset to talk to your friends and stream game audio to your processors.
- Stream game audio to your processors but allow others in the same room to also hear the game audio.





How to stream game audio only to your processors

1. You will need the audio jack cable that is included in the Microphone box.



2. Place one end of the cable into the Microphone.





**SCOTTISH
COCHLEAR IMPLANT
PROGRAMME**

3. If you were using headphones previously, remove these.

4. Plug the other end of the cable either into your controller



Or

into your TV or monitor.






5. Switch the Microphone on.

The indicator light
will turn green
when on.



6. You now need to tell your processor to listen to the Microphone. You can use either the app, remote or processor button to do this.





**SCOTTISH
COCHLEAR IMPLANT
PROGRAMME**

7. You then need to tell the Microphone to listen to the Game audio.

Press the long thin button under the pairing button. You will see a green light flash to the right side. Press the button x3 until the music note lights green.





How to stream game audio to your processor and use a headset to talk to your friends

Follow the steps described on the first few pages, except steps 3 and 4. Do not remove your headset. Plug the other end of the audio cable into the remaining audio slot (either your controller or your TV/monitor). Continue with steps 5 onwards.

8. Go into the console settings. This will differ for each make of console. For Playstation and Xbox you press the console button to access settings.





The following steps are specifically for XBox. The steps for each console will be similar but wording may differ. Speak to a member of the team if you are having difficulties.

9. The menu will appear.
10. Move along to select your 'profile & systems'.
11. Then select settings.
12. Go to General > Volume & audio output.
13. Advanced > Additional options.
14. Party Chat Output > Speakers

How to stream game audio to your processor, use a headset to talk to your friends and allow others to also hear game audio

If you wish for someone else to listen to the game audio also, follow all of the above instructions except for the last step. You will need to plug a headset in also and in Party Chat Output settings select 'Speaker and Headset'. The other user will wear the headphones.